



alpha

OPERATING MANUAL

Issue MD107368
February 2008

Micro Direct, Inc.
803 Webster Street
Lewiston, ME 04240
1-800-588-3381
www.micro-direct.com

Indications for Spirometry

Spirometry has been used extensively to measure lung function capability and to recognize and treat many diseases associated with the impairment of healthy lung functions. Spirometry today provides great insight into the status of any person's health.

Generally speaking, spirometry is a simple diagnostic tool used to define a subject's lung condition. The major indications for spirometry are:

- ✓ Dyspnea (shortness of breath)
- ✓ Exercise induced coughing
- ✓ Chest tightness
- ✓ Smokers over 45 years of age (NLHEP recommendations)
- ✓ Obesity
- ✓ Pre-operative testing
- ✓ Occupational exposure to dust and/or chemicals
- ✓ Ongoing assessment of patients receiving bronchodilator treatments
- ✓ Determination and/or documentation of pulmonary disability
- ✓ Asthma diagnosis
- ✓ Pre-existing pulmonary disease
- ✓ Frequent colds
- ✓ Assessment of congestive heart failure

CPT Codes for Spirometry

94010 - Spirometry Complete

Includes graphic record total and timed vital capacity, expiratory flow rate measurement (s) with or without maximal voluntary ventilation

94040 - Bronchodilation Responsiveness

Spirometry as in 94010, pre and post bronchodilator or exercise

94070 - Bronchospasm Provocation Evaluation

Multiple spirometric determinations after bronchodilator with spirometry as in 94010

94150 - Vital Capacity

Total (separate procedure)

94200 - Maximal Voluntary Ventilation

Maximum breath capacity

94375 - Flow Volume Loop

Respiratory Flow Volume Loop

95070 - Inhalation Bronchial Challenge Testing

(Not including necessary pulmonary function tests), with histamine, methacholine or similar compounds.

94464 - Bronchodilator Administration

Demonstration and/or evaluation of patient utilization of an aerosol generator, nebulizer and meter dose inhaler or IPPB device

Table of Contents

DESCRIPTION OF THE MICRO DIRECT ALPHA	1
MAIN COMPONENTS OF THE MICRO DIRECT ALPHA	1
FEATURES OF THE MICRO DIRECT ALPHA	1
GETTING THE MICRO DIRECT ALPHA READY FOR USE	2
POWER MANAGEMENT IN THE MICRO DIRECT ALPHA	2
BATTERY PACK	2
BATTERY LOW DETECT	3
OPERATING THE MICRO DIRECT ALPHA	3
ENTERING SUBJECT INFORMATION	3
PERFORMING A TEST SESSION	4
<i>Checks to Make before Performing a Test Session</i>	4
<i>Performing a VC Test</i>	5
<i>Performing an FVC Test</i>	7
SAVING THE TEST SESSION.	8
PERFORMING A POST TEST SESSION	8
PRINTING THE TEST SESSION	9
CLEAR RESULTS	9
CHECKING ACCURACY	10
<i>When to Check Accuracy</i>	11
VIEWING CURRENT SUBJECT DETAILS	11
DELETING STORED SUBJECTS/TEST RESULTS	11
CONFIGURATION OPTIONS	12
<i>Preferences</i>	12
<i>Date/Time</i>	13
<i>Report Settings</i>	13
<i>Syringe/Calibration</i>	14
<i>Sound Volume</i>	16
FITTING A NEW PAPER ROLL	16
CLEANING INSTRUCTIONS	17
CLEANING AND DISINFECTING THE MICRO DIRECT ALPHA	17
DISASSEMBLING THE FLEISCH FLOWHEAD	19
REASSEMBLING THE FLEISCH FLOWHEAD	20
FAULT FINDING GUIDE	21

CUSTOMER SERVICE	23
CONSUMABLES AND ACCESSORIES	23
EXPLANATION OF SYMBOLS	24
OTHER LABELS	24
TECHNICAL SPECIFICATIONS	24
CE NOTICE	25
FDA NOTICE	29
DECLARATION OF CONFORMITY	30
GUARANTEE	31




Figure 1

DESCRIPTION OF THE MICRO DIRECT ALPHA

The Micro Direct ALPHA is a desktop spirometer designed for use by trained professionals in the doctor's office, clinic, hospital department, etc. for measuring and archiving tests on human subjects. Demographic data are uploaded or entered via a keypad and stored, together with spirometry test data. Current test data can be viewed on the LCD and printed and downloaded to a PC. There are a variety of backup and other configuration options.

Information about the software can be obtained from the About box. This information can be used if any queries are made to Micro Direct or a service agent.

To access the About box:

1. Press 8 (Configuration) from the Main Menu.
2. Press 6 (About the ALPHA)

MAIN COMPONENTS OF THE MICRO DIRECT ALPHA

- A Micro Direct ALPHA device
- B LCD
- C Keypad
- D Flowhead
- E Printer

FEATURES OF THE MICRO DIRECT ALPHA

The Micro Direct **alpha** features include:

- Very high accuracy, linearity and stability
- Extremely simple to operate keypad with extra-large buttons
- Integral printer for instant results or printing to the Micro Direct Report Utility
- Storage of tests and demographic information
- Ability to exchange information with a PC
- Robust, proven Fleisch-type flowhead. No moving parts or delicate sensors
- Simple cleaning. The flowhead completely dismantles for washing, disinfecting or even autoclaving.

GETTING THE MICRO DIRECT ALPHA READY FOR USE

1. Attach the flowhead to the Micro Direct ALPHA by the dual silicone tubing (Flowhead Connection Tube). Ensure that the colored/ribbed tapping on the flowhead is connected to the ribbed side of the connection in the housing.
2. Open the printer door and feed the paper through the printer (See section Fitting a New Paper Roll).
3. Only use the Micro Direct ALPHA with the purpose-built low voltage PowerSAFE power supply unit with which it is supplied. Attempted use with other power sources may cause irreparable damage and invalidate the warranty. The output from the PowerSAFE is 12 volts DC.
4. Connect the jack plug from the PowerSAFE into the socket on the rear of the Micro Direct ALPHA. Plug the mains plug into a suitable socket, operate the On/Off switch on the rear of the instrument and the Micro Direct ALPHA is ready for use.
5. For portability the Micro Direct ALPHA comes fitted with rechargeable batteries, which allows the device to be used for a period of time without the mains connected.

POWER MANAGEMENT IN THE MICRO DIRECT ALPHA

The Micro Direct ALPHA can be powered using the purpose-built low voltage PowerSAFE unit with which it is supplied or from the internal Battery Pack. When powered from the low voltage PowerSAFE the LED on the front face on the device will be green. The LED will be orange when the device is powered from the Battery Pack.

Battery Pack



The Micro Direct ALPHA is fitted with a rechargeable Battery Pack. This allows the device to be used without the 12V PowerSAFE connected. The battery pack can be re-charged by plugging in the 12 V PowerSAFE. To fully re-charge switch off the Micro Direct ALPHA and leave it plugged in overnight.

The battery pack can also be re-charged by connecting the device to a PC via the USB cable. The USB connector is located on the right side of the device. Switch off the Micro Direct ALPHA when re-charging the batteries via USB.

Note: Operating the Micro Direct ALPHA device from the USB power is not possible. Output from different USB ports vary by brand of PC.

Battery Low Detect

The Micro Direct ALPHA has a number of battery power messages:

	<p>When the Battery Pack starts to run low the <i>Battery Low</i> icon will flash on and off on the Main Menu screen. You will be allowed to continue to use the device. It is advised that you plug in the 12V PowerSAFE to re-charge the batteries and continue testing.</p>
	<p>When the Battery Pack is approaching fully discharged the <i>Battery Discharged</i> icon will appear continuously on the Main Menu screen. If you press any of the keys a warning message will appear and you will not be allowed to proceed into any other screen. Plug in the 12V PowerSAFE to re-charge the batteries and continue testing.</p>

OPERATING THE MICRO DIRECT ALPHA

Entering Subject Information

After turning on the device, you are presented with the **New Subject** screen. A series of dialog boxes appear in which you enter information using the keypad. Instructions appear at the bottom of the screen to guide you through the information to be entered. Once you are ready to move onto the next dialog box, press the 'Enter' key on the keypad. These dialog boxes appear one at a time in the following order:

1. Reference Number
2. Age
3. Height
4. Gender
5. Population Group (on some models)

When entering subject information, you can specify the population group of the subject. There are a number of population groups to select from. These population groups have an associated predicted normal set and correction factor. These normal sets are usually based on sex, age and height and are

useful when comparing a subject's test results with predicted normal values from a suitable reference population.

6. Weight (optional)
7. Smoking History (optional)

If you wish to change the data entered press the 'Delete' key and enter the data again. Pressing the 'Delete' key repeatedly will step back through the dialogue boxes. Once you have finished entering the data press the 'Enter' key to bring you to the **Main Menu**.

The Population Group, Weight and Smoking History fields may not appear in the **New Subject** screen. This is dependent on the variant of the device.

Note: The Subject Information will now be printed to the internal printer if it is selected.

If you wish to enter different subject details select the 'New Subject' option from the **Main Menu** using the keypad. You will be asked if you wish to clear results.

- Press Yes to bring you to the **New Subject** screen, and enter the New Subject details as outlined above.
- Press No to bring you back to the **Main Menu**.

Performing Test Sessions Before Performing a VC Test Session

Note: You can view the results as either a Volume/time graph or a Volume bar graph:

Press key 1 for Volume/time (V/t) graph.

Press key 2 for Volume (V) bar graph.

The VC test can be performed using two methods. Give and demonstrate instructions to the subject so that testing is performed properly:

Method 1. Single Breath VC (recommended):

- a. Inhale as deeply as possible with the flowhead well away from the mouth, insert the mouthpiece, clamping it between the teeth and closing the lips around it to create a perfect seal.
- b. Exhale normally for as long as possible.

Method 2. Multi Breath IVC:

Note: Use the Volume/Time (V/t) graph as the display option for this test method. Also using a SpiroSafe filter or regular cardboard mouthpiece is essential.

- a. Insert the mouthpiece, clamping it between the teeth and closing the lips around it to create a perfect seal.
- b. Breathe in and out normally through the flowhead. This is tidal breathing.
- c. When happy that the subject has achieved steady tidal breathing, continue
- d. Exhale as much as possible
- e. Inhale as much as possible (speed is not important) and when fully inhaled.
- f. Return to tidal breathing, i.e. breathing in and out normally again.
- g. End test by removing mouthpiece.

Before starting a VC test session:

Ensure that the test subject is fully prepared: posture (sitting is recommended); removable false teeth have been removed; tight clothing is loosened and proper instruction on the test is given - including a demonstration by the operator. Fit a new SpiroSafe filter, MicroCheck mouthpiece or regular cardboard mouthpiece and finally fit a disposable nose clip.

Performing a VC Test

Perform the VC test as follows:

1. Select the 'VC Test' option from the **Main Menu** using the keypad.
2. Wait for the 'Exhale to Begin' icon to appear.



This indicates that the Micro Direct ALPHA unit is ready to accept a blow.

The VC values recorded for the blow are tabulated. The best VC value for the current session and the Lower Limit of Normality (LLN) are also displayed.

2. Repeat to perform another VC test if required.

3. After performing the VC tests press the 'Enter key' to exit the **VC Test** screen. This brings you back to the **Main Menu**.

Before Performing an FVC Test Session

Note: You can view the results as either a Volume/time (V/t) or a Flow/Volume (F/V) graph:

Press key 1 for Volume/time graph.

Press key 2 for Flow/Volume graph.

The FVC test can be performed using two methods. Give and demonstrate instructions to the subject so that testing is performed properly:

Method 1. Single Breath FVC (recommended):

- a. Inhale as deeply as possible with the flowhead well away from the mouth, insert the mouthpiece, clamping it between the teeth and closing the lips around it to create a perfect seal.
- b. Exhale as hard and fast as possible for as long as possible.
- c. Optionally at the end of expiration, inhale fully.

Note: It is recommended to use the Volume/Time (V/t) graph as the display option for tests involving inspiration. Using a SpiroSafe Filter or regular cardboard mouthpiece is also essential.

Method 2. Multi Breath FVC:

Using a SpiroSafe filter or regular cardboard mouthpiece is essential for multi-breath.

- h. Insert the mouthpiece, clamping it between the teeth and closing the lips around it to create a perfect seal.
- i. Breathe in and out normally through the flowhead. This is tidal breathing.
- j. When happy that the subject has achieved steady tidal breathing, continue
- k. Exhale as much as possible
- l. Inhale as much as possible (speed is not important) and when fully inhaled
- m. Exhale as hard and fast as possible for as long as possible.
- n. Optionally, return to tidal breathing.

Before starting an FVC test session:

Ensure that the test subject is fully prepared: posture (sitting is recommended); removable false teeth have been removed; tight clothing is loosened and proper instruction on the test is given - including a demonstration by the operator. Fit a new SpiroSafe filter, MicroCheck mouthpiece or regular cardboard mouthpiece and finally fit a disposable nose clip.

Performing an FVC Test

1. Select the 'FVC Test' option from the **Main Menu** using the keypad.
2. Wait for the 'Exhale to Begin' icon to appear.



This indicates that the Micro Direct ALPHA unit is ready to accept a blow.

3. Perform FVC test.
4. The FVC, FEV1, FEV1ratio and PEF values recorded for the blow are tabulated. The best FVC, FEV1, FEV1ratio and PEF for the current session are displayed.

The test quality (QA) is shown at the bottom of the test screen.

The number of tests performed and the Test Grade are shown in the V/t screen. Each test series is graded in relation to its repeatability between acceptable maneuvers. The quality Grades are A, B, C, D and F.

The repeatability (Within) of FVC and FEV1 are shown in the F/V screen. The repeatability information is displayed if at least two tests are performed. I Bars on the F/V graph are shown for FEF25, FEF 50 and FEF75. An I Bar for FVC is also shown on the Volume axis. The upper mark on the I Bars indicates the predicted value for the subject. The lower mark on the I Bar indicates the LLN value for the subject. The I Bars are based on the predicted sets and will be shown if sufficient subject demographics information is entered.

5. Repeat to perform another FVC Test if required – 3 good ‘blow attempts’ are recommended.
6. After performing the FVC tests press the ‘Enter’ key to exit the **FVC Test** screen. This brings you back to the **Main Menu**.

Printing the Test Session

Print the current test session by selecting ‘Print’ from the **Main Menu**.

--oo00oo--

Performing a Post Test Session

A Post test session can be performed on an FVC test session following the administration of a bronchodilator or a bronchial provocation, such as exercise. Post test performance is measured versus pre.

To perform a Post test:

1. Select ‘Post Mode’ from the **Main Menu** using the keypad.
2. If you want to perform a Post Test on the Pre-test Session just performed select ‘Set Post Mode for Current Subject’. This will return you to the **Main Menu** screen. The text Post Mode will appear on the top right hand corner of the screen.
3. If you wish to perform a Post Test on a different subject or Pre-test:
 - a. Select ‘Select a Subject from the Pre-Store’.
 - b. A message ‘Warning! Current Results will be cleared’ will appear. Select Yes and the **Select Subject** screen will appear.
 - c. Select the number of the Pre-test you wish to perform the Post test on. This will return you to the **Post Mode** screen.
 - d. Press ‘Enter’ to return you to the **Main Menu** screen.
4. Perform the Post FVC test as outlined in section Performing a Test Session.

Note: If there is no Pre-Test performed, the message No FVC Tests Performed will appear.

Temporarily storing the Test Session for a later Post Test.

1. Select ‘Post Mode’ from the **Main Menu** using the keypad.
2. Select ‘Save Current Subject Pre Test’ from the **Post Mode** screen using the keypad. A message will appear informing you of the

memory location where the test session will be saved. There are nine memory locations on the Micro Direct ALPHA.

Note: If a suitable Compact Flash card (available as an accessory from Micro Direct) is inserted into the Compact Flash connector at the right side of the device, then all test blows and not just the best three will be saved to the Compact Flash card. Results are saved as per the format outlined in the European Respiratory Journal, 2005; 26: Pages 319-338: ATS/ERS Task Force: Standardisation of Lung Function Testing.

About Printing the Test Session

You can print the current test session for the subject by selecting 'Print' from the **Main Menu**.

The Micro Direct ALPHA has an internal printer but it can also be connected through the USB port at the side of the unit to the Micro Direct Reports Utility, so that the report can be written to a PC or a printer on the network.

The information printed on the session report can be configured to suit individual requirements. Refer to section on [Report Settings](#).

The test parameters on the report will vary according to regional requirements. The test parameters and their definitions are available from the Help files, which can be accessed from the **Main Menu**.

Note: The internal (thermal) printout will fade over time when exposed to light or heat. If a permanent record is required, photocopy or scan the thermal printout or still better send the report to the Micro Direct Reports Utility. This creates a PDF image of the spirometry report, which can be printed, re-named, e-mailed, etc. Refer to the section on [Configuration](#) for information on selecting the internal or Micro Direct Reports Option.

Clearing Results

If you wish to delete the current session you can do this as follows:

1. Select the 'Clear Results' option from the **Main Menu** using the Keypad.
2. A message will appear 'Warning! Current Results will be cleared. Do you wish to proceed?' in the **Clear Results** screen. Select

'Yes' to delete the current results and return to the **Main Menu**.
Select 'No' to cancel the delete and return to the **Main Menu**.

Checking Accuracy

All spirometry standards (e.g. ATS/ERS/BTS/ANZRS) recommend checking the accuracy of lung function measuring devices at least daily with a 3-L syringe to validate that the instrument is measuring accurately. The Micro Direct ALPHA should never be outside accuracy limits unless damaged or in a fault condition. In this event, see the [fault-finding guide](#). In normal use, calibration traceability certification is recommended as a part of the routine annual service.

ATS/ERS recommendations require that the difference between the volume measured by the spirometer and the volume pumped into the spirometer from a syringe is within 3%.

Follow these steps to check the accuracy of the unit.

1. Select Accuracy Check from the **Main Menu** using the keypad.
2. Pump air through the flowhead to bring it to ambient temperature.
If the flowhead has very recently been used for testing or has come from a cold environment, its temperature should be equilibrated with ambient by pumping air through it from the syringe several times.
3. Press the 'Enter' key to bring you into the **Accuracy Check** screen and follow the on-screen instructions.

*Note: Press the 'Del' key to exit the **Accuracy Check** screen and return to the **Main Menu**. The accuracy check will not be logged to the Micro Direct ALPHA memory in this case.*

4. If an Accuracy Check report is required select the Report option.

Note: If the device is outside calibration you will be given the option to update the calibration. If you select this option you will be brought through the accuracy check routine again. To adjust the calibration it is necessary to enter the 'configuration' option.

When to Check Accuracy

- In accordance with your own established procedures
- After annual maintenance checks
- After cleaning or disassembling spirometer for any reason
- After adjusting calibration
- If the flowhead or tubing has been dropped or damage is suspected.

Viewing Current Subject Details

To view the Current Subject details:

1. Select 'Post Mode' from the **Main Menu** using the keypad.
2. Select 'View Current Subject' from the **Post Mode** screen.
3. The **Current Subject** screen will be displayed.

Note: You will not be able to modify any of the subject details.

Deleting Stored Pre-Test Subjects/Test Results

To delete individual Stored Subjects/Test results:

1. Select 'Post Mode' from the **Main Menu** using the keypad.
2. Select 'Delete Stored Tests' from the **Post Mode** screen.
You are presented with the list of the Subjects and the associated test stored on the device. A maximum of nine Subject/Tests can be stored on the device in locations 1-9 as indicated in the **Delete Subject** screen.
3. Select the number of the Subject/Test you wish to delete using the keypad. This location will then be marked as 'Empty'.
4. Continue to delete test as required.
5. Press 'Enter' to return to the **Post Mode** screen.

To delete all Stored Subjects/Test results:

1. Select 'Post Mode' from the **Main Menu** using the keypad.
2. Select 'Clear Pre Store' from the **Post Mode** screen.
3. A message 'Warning! All data in pre store will be cleared. Do you wish to proceed?' will appear in the **Clear Results** screen.
4. Select 'Yes' to delete the all data and return to the **Post Mode** screen. Select 'No' to cancel the delete and return to the **Post Mode** screen.

Configuration Options

There are a number of Configuration options available on the Micro Direct ALPHA device. To access these, press the 'Configuration' option on the **Main Menu** using the keypad. The options available are:

Preferences

This will allow you to configure the device to your own requirements. Selecting this will give you the following options:

Printer

This allows you to either print the test report to the Internal Printer or send it to the Micro Direct Reports Utility on a PC. Press key '1' on the keypad to switch between the two options.

Note: In order to send the report to the Micro Direct Reports Utility it is necessary to have the Micro Direct Reports Utility installed on your PC and the Micro Direct ALPHA connected to your PC via a USB cable.

Units

This allows you to select the units as Metric or US (Imperial). Press key '2' on the keypad to switch between the two options.

Graph Type

This allows you to select the graph type to be shown as default in the FVC test screen, Press key '3' on the keypad to switch between the Volume/Time (V/T) and Flow/Volume (F/V) graphs.

Test Acceptability

This allows you to manually accept the tests performed, or allow the device to determine test acceptability (automatic). Press key '4' to switch between Manual and Automatic.

VC Test Type

This allows you to select the graph type to be shown as default in the VC test screen, Press key '5' on the keypad to switch between the Volume/Time (V/T) and Volume Bar graphs.

Date/Time

1. Select 'Date/Time' from the **Configuration** screen using the keypad.
2. In the **Date/Time** screen press key '1' to change the year. Enter the revised year as required using the keypad and press 'Enter'. The Month, Day, Hour and minute fields are modified in the same way.
3. To modify the Date Format press key '6' to switch between the different date formats:
DD/MM/YYYY
MM/DD/YYYY
YYYY/MM/DD
4. Press 'Enter' to save the change.
5. To modify the Time Format press key '7' to switch between 24 hour and 12 hour.
6. Press 'Enter' to save the change.
7. Once you have finished modifying the date and time settings press 'Enter' to save the changes and return to the **Configuration** screen. Press 'Delete' to cancel the changes made and return to the **Configuration** screen.

Report Settings

The information printed on the session reports can be configured to suit individual requirements. The Micro Direct ALPHA can print to the internal printer or send the report to the Micro Direct Reports Utility.

1. Select 'Report Settings' from the **Configuration** screen using the keypad. You are presented with the list of items that can be configured on the test reports.
2. In the **Report Settings** screen press key '1' to switch the Volume/time (V/t) graph on or off.
3. Press key '2' to switch the Flow/Volume (F/V) graph on or off.
4. Press Key '3' to switch the Comments on or off.
5. Press key '4' to switch the Test QA messages on or off.
6. Press key '5' to switch the Interpretation of the test results on or off.
7. Press key '6' to switch between showing SDS (Standard Deviation Score) or %Predicted in the session table of results. In the session table of results the % of Predicted value will be printed by default. If SDS is on, the SDS (Standard Deviation Score) will be printed instead.
8. Press key '7' to either show the Best Test only on the report or the Three Tests saved to the database.

9. Press key '8' to switch color on or off. This option is not available for the internal printer.
10. When you have finished modifying the Report Settings press 'Enter' to save the changes and return to the **Configuration** screen. Press 'Delete' to cancel the changes made and return to the **Configuration** screen.

Syringe/Calibration

The Micro Direct ALPHA should never be outside accuracy limits unless damaged or in a fault condition. In this event, see the [fault-finding guide](#). In normal use, calibration traceability certification is recommended as a part of the routine annual service.

1. Select the Syringe/Calibration menu using the keypad. You are presented with three options:
 - a. Set Syringe Volume
 - b. Linearity Check
 - c. Calibration

Set Syringe Volume:

- i. Select Precision Syringe from the **Syringe/Calibration** screen using the keypad.
- ii. Follow on screen instructions to enter the volume of the calibrated syringe you are using.
- iii. Press 'Enter' to save the new volume entered and return to the **Syringe/Calibration** screen. Press Delete to cancel the changes made and return to the **Syringe/Calibration** screen.

Linearity Check

- i. Select Linearity Check from the **Syringe/Calibration** screen using the keypad.
- ii. Pump air through the flowhead to bring it to ambient temperature. If the flowhead has very recently been used for testing or has come from a cold environment, its temperature should be equilibrated with ambient by pumping air.
- iii. Press the 'Enter' key to bring you into the **Linearity Check** screen.
- iv. Using a 3L Calibrated syringe pump air into the flowhead at a slow rate of <2L/S. Immediately withdraw the syringe at a slow rate. This

maneuver should show on the graph between the two red lines. If it is a correct maneuver the table on the screen will show 'Test 1', and the FVC and FIVC values will be updated.

*Note: Press the 'Del' key to exit the **Linearity Check** screen.*

- v. Repeat for the slow rate three times in total.
- vi. Repeat the procedure outlined in iv & v for a medium rate >2L/s and <6L/s. This maneuver should show on the graph between the red and green lines. If it is a correct maneuver the test number and the FVC and FIVC values will be updated in the table.
- vii. Repeat for the medium rate three times in total.
- viii. Repeat the procedure outlined in iv & v for a fast rate >6L/s. This maneuver should show on the graph between outside green lines. If it is a correct maneuver the test number and the FVC and FIVC values will be updated in the table.
- ix. Repeat for the medium rate three times in total.
- x. When all the manoeuvres are complete press 'Enter' for the result.
- xi. If a Linearity Check report is required select the Report option.

Calibration

- i. Select Calibration from the **Syringe/Calibration** screen using the keypad.
- ii. Pump air through the flowhead to bring it to ambient temperature. If the flowhead has very recently been used for testing or has come from a cold environment, its temperature should be equilibrated with ambient by pumping air through it from the syringe several times.
- iii. Press the 'Enter' key to bring you into the **Calibration** screen and follow the on-screen instructions.

*Note: Press the 'Del' key to exit the **Calibration** screen.*

- iv. If a Calibration report is required select the Report option.

Sound Volume

1. Select the 'Sound Volume' option from the **Configuration** screen using the keypad.
2. Adjust the sound level up by pressing key '2'.
3. Adjust the sound level down by pressing key '8'.
4. When you have finished changing the sound level press 'Enter' to save the changes and return to the **Configuration** screen.

Fitting a New Paper Roll

The Micro Direct ALPHA is supplied with a roll of paper fitted into the printer. *(Note: The outside of the roll will not give a strong print impression due to loss of sensitivity on exposure to light.)*

To replace the paper:

1. Open the printer door using the finger catch on each side of the door. This will expose the head of the printer and the paper roller mechanism.
2. Take out the empty roll by unclipping the roll support rod from the holder.
3. Put the roll support rod into the new roll of paper and unroll about 15 cm (6 inches) of paper. *(Note: To make fitting easier create a point in the middle of the paper by tearing off the two corners of the leading edge of the paper.)*
4. Lay the paper roll between the back of the Micro Direct ALPHA unit and the inside of the open door, with the paper coming out from the bottom of the roll and pointing toward the printer. This will allow you more access to the printer feed mechanism.
5. In the home screen, lift the green lever on the printer.
6. Feed the leading edge of the paper into the bottom slot of the printer until the paper appears through the top of the printer. The paper can now be pulled through. To aid in feeding the paper through the printer, press the 'Enter' key on the **Main Menu** screen. This will feed a short length of paper through the printer.
7. Close the green lever.

8. Fit the new paper roll onto the holding clips.
9. Hold the paper over the paper tear bar and close the door.

Note: The Micro Direct logo should appear on the right edge of the paper and be facing you.

Warning: Paper tear bar contains sharp edges. Users should take care not to cut/scrap their fingers.

Cleaning Instructions

Cleaning and Disinfecting the Micro Direct ALPHA

A new mouthpiece (SpiroSafe, MicroCheck or regular cardboard mouthpiece) should be used for each subject. A delay of at least 5 minutes should be allowed between subjects to allow settling of previously aerosolized particles in the measuring device.

It is recommended that the flowhead be regularly cleaned according to the guidelines of the user's facility.

In the event of visible contamination of the flowhead cones or flowhead element, they should be cleaned or disinfected as described in the accompanying table. The flowhead conditioning meshes should also be replaced in the event of damage, or if visibly contaminated.

The frequency of cleaning and disinfecting is dependent on the Facility's Risk Assessment, usage, and test environment, but it should be at least monthly or every 100 subjects (500 blows).

It is recommended that the flowhead—flowhead complete and flowhead connection tube—be replaced annually.

Refer to [Figure 2](#) for identification of parts.

Table of Materials Used & Cleaning/Disinfection Methods

This listing of materials used is given to provide users with information to allow the assessment of other cleaning and disinfecting procedures available in the facility on this device.

Part	Material	Clean/Disinfect	Autoclave Possible?	Recommended Disinfectants
Case Exterior	ABS	Clean	No	Wiping with a 70% isopropyl alcohol impregnated cloth provides a suitable form of cleaning and low-level disinfection. This may be preceded by cleaning with an anti-static foam cleaner if necessary. Note: Ensure isopropyl alcohol does not come in contact with the screen. <i>Warning: Paper tear bar contains sharp edges. Users should take care not to cut/scrap their fingers.</i>
White Flowhead Tube	Silicone Rubber	Clean	Viable	
Screen	Electrode with Anti-Newton-Ring Treatment	Clean	No	Lightly wipe the surface with cotton pad or other soft material. NOTE: DO NOT use chemicals such as acetone, toluene, ethanol or isopropyl alcohol. DO NOT wipe in a circular motion. Strokes should be either up/down or over/back.
Fleisch Element	Aluminum, Stainless Steel	Clean	Viable	Disinfect by immersion in sodium dichloroisocyanurate solution at 1000 ppm concentration of free chlorine for 15 minutes (see following section for recommended cleaning/disinfection method for the Micro Direct ALPHA Flowhead) The flowhead may also be disinfected by autoclaving at 134°C for 3 minutes or 120°C for 20 minutes.
Flowhead Body	Aluminum & Acetyl	Clean & Disinfect	No	
Flowhead Cone	TPX	Clean & Disinfect	Viable	
Flowhead End Cap	TPX	Clean & Disinfect	Viable	
Flow Conditioning Meshes	Acetyl and Polyester	Dispose	No	

All external parts of the Micro Direct ALPHA require **cleaning**, i.e. the removal of visible particulate contamination. The parts of the Micro Direct ALPHA that make up the flowhead, which comes into contact


with subjects being tested, also require **disinfecting**. A spirometer is not designed as a 'sterile' device.

Definitions of cleaning and disinfection are as defined in "Sterilization, Disinfection and Cleaning of Medical Equipment: Guidance on Decontamination from the Microbiology Committee to Department of Health Medical Devices Directorate, 1996".

Recommendations for chemical disinfectants are derived from the PHLS publication "Chemical Disinfection In Hospitals 1993".

Disassembling the Fleisch Flowhead

1. Remove the cone and the end cap from the flowhead.
2. Remove the flow conditioning meshes from inside the cone and the end cap, and discard them.
3. To remove the flowhead body from the Fleisch element, place the Fleisch element on a hard, flat surface with the largest diameter at the top. Push down on the flowhead body with thumbs and forefingers until it reaches the flat surface. A final pulling and twisting action will separate the parts.



1. Flowhead Complete – 61030
2. 'O' rings - 2120013
3. Flowhead End Cap -62006SPR
4. Flow Conditioning Meshes - 42084
5. Flowhead Cone - 62019SPR
6. Fleisch Element - 62055SPR
7. Flowhead Body – 61020
8. Lubrication: Silicone Grease – 30961SPR

Figure 2: Flowhead Assembly

4. Clean each separate part of the flowhead by washing in a mild detergent and removing particulate contamination. To clean the Fleisch element, swill vigorously in water with mild detergent or use an ultrasonic bath. Do not attempt to "rub" or "scrub" at capillaries.

The flowhead body (7) does not require disinfection, but may be cleaned/disinfected with the rest of the flowhead for convenience.

5. Rinse all parts in clean water.
6. Disinfect by immersion in sodium dichloroisocyanurate (NaDCC) solution at 1,000 ppm concentration of free chlorine for 15 minutes. Prepare disinfectant solution as directed in the manufacturer's guidelines.
7. Rinse with very hot water to aid later drying.
8. Leave to dry completely before reassembling. Drying the Fleisch element components may require placing them in a warm place overnight. A drying cabinet is ideal.

Always follow the safety guidelines given by the manufacturer of cleaning and disinfectant chemicals.

Reassembling the Fleisch Flowhead

1. Examine the Fleisch element to ensure that no liquid or particles remain in the holes, grooves or pressure tappings.
2. Check the 'O' rings for damage and ensure that they are correctly positioned within the grooves.
3. Apply a very small amount of silicone grease to 'O' Rings and inside the surfaces of the flowhead body. Wipe off any visible amounts of grease. Ensure that the tiny annular holes on the outside of the Fleisch element are not blocked.
4. When re-assembling the flowhead, ensure that the blue pressure tapping is nearest to the largest diameter of the Fleisch element.
5. Ensure that the flowhead body is pushed fully home and rotate it so that the pressure ports are approximately 180° opposite the end of the Fleisch element coil.
6. Fit new flow conditioning meshes to both the flowhead cone and the flowhead end cap.
7. Push the flowhead end cap onto the larger diameter of the Fleisch element and push the flowhead cone onto the smaller diameter.

8. When attaching the flowhead connection tube ensure that the matching colored/serrated edge pressure tapings on the flowhead and the Micro Direct ALPHA are connected to each other.
9. It is recommended that an accuracy check is carried out following reassembly to verify correct operation and accuracy.

FAULT FINDING GUIDE

Problem Fault Symptoms:	<ul style="list-style-type: none"> • Accuracy check variations > +/-3% • False readings suspected
Possible Causes: (In probable order)	<ul style="list-style-type: none"> • Recheck Calibration with reference to section Checking Accuracy • Was the correct syringe volume selected? • An accuracy check is required after cleaning/disinfecting the flowhead Fleisch element assembly. • Flowhead cone Fleisch element filter mesh missing or blocked. • Flowhead body pressure port holes blocked. • Flowhead Fleisch element assembly sealing 'O' rings damaged. • Flowhead Fleisch element assembly not dried thoroughly. • Flowhead Fleisch element assembly blocked. • Flowhead body tubing from pressure ports to main PCB blocked – contact support. • Main PCB failure – contact support.
Problem Fault Symptoms:	<ul style="list-style-type: none"> • Test begins automatically • Volume accumulates automatically without the subject blowing. • Very small VC or FVC test displayed
Possible Causes: (In probable order)	<ul style="list-style-type: none"> • Flowhead and/or tubing not stationary at the start of test. Hold them steady until the 'Ready to Blow' prompt appears. • Return to Main Menu and re-enter the test routine.

Problem Fault Symptoms:	<ul style="list-style-type: none"> •Rocking device.
Possible Causes: (In probable order)	<ul style="list-style-type: none"> •Check for damaged or missing rubber feet. •If any of the rubber feet are damaged or missing replace all six rubber feet.
Problem Fault Symptoms:	<ul style="list-style-type: none"> •Reversed or no volume measurements.
Possible Causes: (In probable order)	<ul style="list-style-type: none"> •Ensure tubing is connected correctly. Ribbed side of the tubing should be connected to the ribbed half of the connector on the Micro Direct ALPHA device and the blue tapping on the flowhead connector •Ensure that the flowhead connecting tube is not pinched or trapped.
Problem Fault Symptoms:	<ul style="list-style-type: none"> •Cannot print to internal printer.
Possible Causes: (In probable order)	<ul style="list-style-type: none"> •Check that internal printer is selected in the Configuration screen. •Check the paper is loaded correctly and not reversed. •Ensure the green flap on the printer is pressed down. •Internal printer failure – contact support.
Problem Fault Symptoms:	<ul style="list-style-type: none"> •Cannot print to PC (Micro Direct Reports Utility). •Corrupt or missing data on printout.
Possible Causes: (In probable order)	<ul style="list-style-type: none"> •Check that external printer is selected in the Configuration screen. •Check USB cable is connected between Micro Direct ALPHA and the PC. •Check to ensure the Micro Direct Reports Utility is correctly installed. •Check to ensure the required software drivers are installed on the PC. •Main PCB failure – contact support.
Problem Fault Symptoms:	<ul style="list-style-type: none"> •Cannot read screen.

Possible Causes: (In probable order)	<ul style="list-style-type: none"> •Ensure the switch on the back of the unit is in the 'On' position. •LCD failure – contact support. •Main PCB failure – contact support.
---	--

CUSTOMER SERVICE


Service and repairs should be carried out only by the manufacturer, the approved importer or by Service Agents specifically approved by Micro Direct.

For the names and addresses of approved Micro Direct Service Agents or to arrange for an inservice via the telephone, please refer to the contact information at the start of this manual.



CONSUMABLES AND ACCESSORIES

Cat. no	Description
3395	MicroCheck Mouthpieces (200)
3304	Disposable Nose clips (25)
3385	SpiroSafe Filters (100)
3325	3-L Precision Syringe
42084	Flow Conditioning Mesh (10)
66149	Thermal Printer Paper (5)
67252	USB Cable
41195	12V DC PowerSAFE
41198	2 Pin Mains Input Module (USA)
61030	Flowhead Complete
42029SPR	Flowhead Connection Tube
65354	CD with User Manual
65049	Test Data Storage Card
65030SPR	Micro Direct Reports Utility

EXPLANATION OF SYMBOLS

k	Type BF equipment
j	Class II
VA	Power rating
V 	Voltage DC
h	Attention (reference relevant section in manual)

OTHER LABELS

	Power input connector
	USB connector

TECHNICAL SPECIFICATIONS

Product	Micro Direct ALPHA
Model	6000
Flow detection principle	Fleisch type pneumotachograph
Back pressure	Less than 0.1kPa/L/second @ 14L/s, complies with ATS/ERS 2005
Volume detection	Flow integration sampling @ 100Hz
Maximum test duration	90 seconds
Maximum displayed volume	10 L
Volume accuracy	Better than $\pm 3\%$
Voltage/Frequency	110-250 V; approximately 50/60 Hz
Accuracy when operated in operating temperature range conditions	Flow $\pm 10\%$ Max. flow rate ± 16 L/s Min. flow rate ± 0.02 L/s
Operating temperature range	ATS/ERS limits: 17–37°C Design limits: 10–40°C

Performance standards the Micro Direct ALPHA meets or exceeds	ATS/ERS 2005 & ISO 23747:2007
Safety standards	EN ISO 60601
QA/GMP standards	EN ISO 13485:2003; FDA 21CFR820; CMDR; JPAL
Size	12" x 10" x 3"
Weight	4 pounds 7 ounces
Parameters measured	Varies by country variant
Printer	Thermal
Communications	USB and Compact Flash


Note: All values displayed by the Micro Direct ALPHA are expressed as BTPS values.

CE NOTICE

Marking by the symbol **I** indicates compliance of the Micro Direct ALPHA to the Medical Devices Directive of the European Community. Such marking is indicative that the Micro Direct ALPHA meets or exceeds the following technical standards:

Guidance and manufacturer's declaration – electromagnetic emissions		
The ALPHA is intended for use in the electromagnetic environment specified below. The customer or the user of ALPHA should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The ALPHA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The ALPHA is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage PowerSAFE network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Class A	

Guidance and manufacturer's declaration – electromagnetic immunity			
The ALPHA is intended for use in the electromagnetic environment specified below. The customer or the user of the ALPHA should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for PowerSAFE lines ±1 kV for input/output lines	±2 kV for PowerSAFE lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on PowerSAFE input lines IEC 61000-4-11	<5 % 100V (>95% dip in 100V) for 0,5 cycle 40 % 100V (60% dip in 100V) for 5 cycles 70 % 100V (30% dip in 100V) for 25 cycles <5% 100V (>95% dip in 100V) for 5 sec	Performance A Performance A Performance A Performance A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ALPHA requires continued operation during power mains interruptions, it is recommended that the ALPHA be powered from an uninterruptible PowerSAFE or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not Applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity			
The ALPHA is intended for use in the electromagnetic environment specified below. The customer or the user of the ALPHA should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80 MHz in ISM bands	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the ALPHA including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3V/m from 80MHz to 2.5 GHz	$d = 1.2\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = 2.3\sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Recommended separation distances between portable and mobile RF communication equipment and the ALPHA			
The ALPHA is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ALPHA can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ALPHA as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.1m	0.1m	0.2m
0.1	0.4m	0.4m	0.7m
1	1.2m	1.2m	2.3m
10	3.7m	3.7m	7.4m
100	11.7m	11.7m	23.3m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Medical Devices may be affected by cellular telephones and other personal or household devices not intended for medical facilities. It is recommended that all equipment used near the Micro Direct product comply with the medical electromagnetic compatibility standard and to check before use that no interference is evident or possible. If interference is suspected or possible, switching off the offending device is the normal solution, as is required in aircraft and medical facilities.


Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided,

Portable and mobile RF communications equipment can affect medical electrical equipment.

FDA NOTICE

Caution: Federal Law restricts this device to sale by, or on the order of a physician.

DECLARATION OF CONFORMITY


Product: 
Alpha



Micro Direct hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

- European Medical Devices Directive {MDD} 93/42/EEC.

This device, classified as 2a as per Annex IX of MDD 93/42/EEC, meets the following provisions of Annex II of the Medical Devices Directive as per Article 11, section 3a, excluding point 4 of Annex II.



- Canadian Medical Device Regulation {CMDR}
- FDA Quality System Regulation {QSR} 21 CFR 820.
- EN ISO 13485: 2003. Medical devices. Quality management systems. Requirements for regulatory purposes.

Certifying Body {for 93/42/EEC and CMDR}: British Standards Institute {BSI}

Certificate Nos. CE 00772, MD 82182, FM 83550

Signed on behalf of Vitalograph (Ireland) Ltd.



B. R. Garbe.
Group Managing Director

GUARANTEE

Subject to the conditions listed below, Micro Direct Inc. and its associated companies, (hereinafter called the Company) guarantee to repair or at its opinion replace any component thereof, which, in the opinion of the Company is faulty or below standard as a result of inferior workmanship or materials.

The conditions of this Guarantee are:

1. This Guarantee shall only apply to hardware defects which are notified to the Company or to its accredited distributor within 1 year of the date of purchase of the equipment, unless otherwise agreed in writing by the Company
2. Software (meaning computer software or user installable modules) is guaranteed for 90 days from the date of purchase.
3. The Company warrants that the software when correctly used in conjunction with the hardware will perform in the manner described in the Company's literature and user manuals. The Company undertakes to rectify at no expense to the customer any software failure notified within the period stated above, provided that the failure can be recreated and the software has been installed and used in accordance with the user manual. Notwithstanding this clause, the software is not warranted to be free of errors.
4. This Guarantee does not cover any faults caused by accident, misuse, neglect, tampering with the equipment, use of consumable items or parts not approved by the Company, or any attempt at adjustment or repair other than by personnel accredited by the Company, nor does it cover reinstatement of any configuration changes caused by the installation of any software.
5. If a defect occurs, please contact the supplier from whom it was purchased for advice. The Company does not authorize any person to create for it any other obligation or liability in connection with Micro Direct equipment
6. This Guarantee is not transferable and no person, firm or company has any authority to vary the terms or conditions of this Guarantee.
7. To the maximum extent permitted by law, the Company does not accept liability for any consequential damages arising out of the use of, or inability to use any Micro Direct equipment.
8. This Guarantee is offered as an additional benefit to the Consumer's statutory rights and does not affect these rights in any way.

ICD-9 Codes for Spirometry

Diagnosis	ICD-9-CM Codes
Smokers over 40	491.0
Shortness of Breath	518.82
Chronic Cough	464.4, 493.9
Frequent Coughs	460 or 465, 465.0, 465.8, 465.9
Allergic Rhinitis	506, 506.0, 506.1, 506.2, 506.3, 506.4, 506.9
Occupational Exposure to Dust or Chemicals	506, 506.0, 506.1, 506.2, 506.3, 506.4, 506.9
Scoliosis	737, 737.0, 737.1, 737.10, 737.12, 737.19, 737.2, 737.20, 737.21, 737.22, 737.29, 737.3, 737.30, 737.31, 737.32, 737.33, 737.34, 737.39, 737.4, 737.40, 737.41, 737.42, 737.43, 737.8, 737.9
Pigeon Chest	738.3, 754.82
Barrel Chest	783.3
Diagnosis of Asthma	493, 493.0, 493.1, 493.2, 493.9
Diagnosis of Bronchitis	491, 491.0, 491.1, 491.2, 491.8, 491.9
Diagnosis of other COPD	496
Pre-Operative Evaluation	518.5
Wheezing	786.09
High Risk Medication	V58.69