



The PRIMACS^{SNC-100} Total Carbon & Total Nitrogen Analyzer

The analysis of Carbon and Nitrogen is important in soil, plant, animal feed, food samples, sediments and sludge samples. The Skalar Primacs^{SNC-100} analyzer provides fast and reliable analysis of Total Nitrogen (TN)/Protein, Total Organic Carbon (TOC), Total Elemental Carbon (TEC) and Total Inorganic Carbon (TIC) in solid and liquid samples.

The Primacs^{SNC-100} combines the analysis of TC and TN in one unit. High temperature combustion with Non Dispersive Infrared detection (NDIR) is used for the analysis of TOC, TEC and TIC. The temperature settings are variable to get optimum combustion for different sample matrices and to allow the analysis of TEC. The determination of TN/Protein is based upon DUMAS methodology and detection with Thermal Conductivity (TCD). TIC can also be analyzed using automatic acidification and purging.




Operating Principle

Combustion

The samples are weighed in re-usable crucibles and placed in the autosampler. In the high temperature combustion furnace, at 1200°C, Carbon is completely oxidized to CO₂ and Nitrogen is converted into N_xO_y in the presence of oxygen. In the oxidation/reduction oven all Nitrogen is reduced to N₂. First CO₂ is measured by NDIR. Secondly the N₂ gas is measured by TCD.

Acidification

The samples are introduced in the IC reactor. Acid is automatically added to the reactor. IC is then released as CO₂ and measured by the NDIR detector.



Autosampler

The Primacs^{SNC-100} contains a large integrated 100 positions autosampler with sample weights up to 3 g of solid material. Due to the unique vertical sample introduction, the sample ashes remain in the crucible after the analysis and are taken out of the instrument with removal of the crucible.

Quality Control

The analyzer uses various control systems to guarantee correct operation and accurate results such as:

- An internal active temperature stabilization system, which eliminates influences of room temperature
- An automated control system, which checks the gas flow and pressure on various places in the system

Software

The software displays the carbon and nitrogen peaks simultaneously in real-time and the results can easily be printed or exported to a LIMS system. Whenever priority samples have to be analyzed, the work list can be extended during the run.

The Primacs^{SNC-100} complies with international regulations such as - CEN, ISO 10694, NEN-EN 131137, AOAC 990.03, AOAC 992.15, AACC 46-30 and ASBC.

General Characteristics

Analytes	Total Carbon - Total Nitrogen-Protein, TIC, TEC and TOC
Analysis Method	High temperature catalytic combustion, TN according to DUMAS method, TIC by Acidification, TIC-TEC-TOC by High temperature combustion with controlled temperature profile
Detection Method	TN-Protein by Thermal Conductivity Detection (TCD) Total Carbon, TIC, TEC and TOC by Non Dispersive Infrared Detection (NDIR)
Applications	Solid and liquid samples in - Food, soil, plant, milk, animal feed, etc. Soil (TIC, TEC, TOC by Combustion with temperature profile)
Auto sampler	Random Access, 100 positions, vertical sample introduction
Sample introduction	Unique vertical "bottom-to-top" sample introduction system
Features	Automatic balance interfacing Automatic leak test Back flush system to remove ambient air for accurate low level analysis Re-usable quartz crucibles Complies with international regulations such as ISO, EN, AOAC, ASBC etc.
Data processing	Area calculation (multi-point linear regression), interfaced balance, raw data storage and connection to LIMS.
Output:	I/O: Windows based SNCAccess software
Results output:	Computer screen, Harddrive, CD-Rom, USB-Drive, printer and ASCII files

Operational and Performance Characteristics

Measuring range	Carbon: 0.01 – 200 mg C Absolute, Nitrogen: 0.02-100 mg N absolute
Min. detection limit:	Nitrogen: 0.02 mg N absolute, Carbon: 0.01 mg C absolute
Analysis time	Approx. 3 - 5 minutes for TC/TN or TIC
Sample size	Solid samples: up to 3 gram (100-1000 mg nominal) Liquid samples: up to 1000 mg (µl) max.
Furnace Temperature	Combustion oven 1200° C 2-Zone Oxidation/Reduction oven 750/600° C TIC reactor at 150° Celsius
Reproducibility	≤ 1 % RSD for Carbon ≤ 0.5 % RSD for Nitrogen

Physical Characteristics

Gas	Carrier gas; Helium gas 99.99 % pure, at 300 kPa. Combustion gas; Oxygen gas 99.99 % pure, at 300 kPa.
Power requirements	220V – 240 V/50/60Hz.
Dimensions (hxdxw)	99 x 49 x 90 cm, approx.
Weight	120 kg, approx.

Skalar's Headquarters
Skalar Analytical B.V.
Tinstraat 12
4823 AA Breda
The Netherlands
Tel. +31 (0)76 5486 486
Fax. +31 (0)76 5486 400



ISO9001 Certified
ISO14001 Certified

USA Headquarters
Skalar, Inc.
5012 Bristol Industrial Way # 107
Buford, GA 30518
Tel. + 1 770 416 6717
Toll Free: 1 800 782 4994
Fax. + 1 770 416 6718

Asia / Middle East
Skalar Analytical India Pvt Ltd
No. 7/4, Pappathiammal Street
Jain Colony, Kodambakkam
Chennai - 600024 - India
Tel. + 9144 2483 7007
Fax. + 9144 2483 6006

Germany
Skalar Analytic GmbH
Gewerbestraße Süd 63
41812 Erkelenz
Germany
Tel. + 49 (0)2431 96190
Fax. + 49 (0)2431 961970

United Kingdom
Skalar (UK) Ltd.
Breda House,
Millfield Industrial Estate,
Wheldrake, York, YO19 6NA
Tel. + 44 (0)1904 444800
Fax. + 44 (0)1904 444820

France
Skalar Analytique S.A.R.L.
35 - 37, rue Berthollet
94110 Arcueil
France
Tel. + 33 (0)1 4665 9700
Fax. + 33 (0)1 4665 9506

Canada
Skalar, Inc.
Unit # 200, 270 Orenda Road
Brampton, L6T 4X6
Toll Free: 1 800 782 4994
Tel. + 1 770 416 6717
Fax. + 1 770 416 6718

Email: info@skalar.com
©Copyright Skalar 2015

Internet: www.skalar.com
Publication No. 0406012A. US

Skalar reserves the right to change specifications and appearance of the equipment without further notification.